

Anti-IKAP antibody (STJ501408)

STJ501408

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-IKAP is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Chicken
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	230233
Gene Symbol	Elp1
Uniprot ID	ELP1_MOUSE
Immunogen	Recombinant IkappaB kinase complex-associated protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-IKAP Antibody (STJ501408) of recombinant IKAP protein. Concentrations of recombinant IKAP protein are listed. Antibody at 1:500 dilution in buffer.