

Anti-IKAP antibody (1280-1330) {FITC} (STJ501410)

STJ501410

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IKAP (1280-1330) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	230233
Gene Symbol	Elp1
Uniprot ID	ELP1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1280-1330 on mouse Elongator Complex Protein 1.
Immunogen Region	1280-1330
Specificity	
Immunogen Sequence	