

## Anti-IL17RA antibody (110-150) {FITC} (STJ501424)

STJ501424

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-IL17RA (110-150) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IF, IHC, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.65-0.95 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">23765</a>
<b>Gene Symbol</b>	<a href="#">IL17RA</a>
<b>Uniprot ID</b>	<a href="#">I17RA_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 110-150 on human Interleukin-17A precursor protein.
<b>Immunogen Region</b>	110-150
<b>Specificity</b>	This antibody is specific for Interleukin-17A and will not cross react with any other Interleukin-17 variants.
<b>Immunogen Sequence</b>	