

Anti-IL18RAP antibody (500-550) {Biotin} (STJ501438)

STJ501438

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IL18RAP (500-550) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8807
Gene Symbol	IL18RAP
Uniprot ID	I18RA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human Interleukin-18 receptor accessory protein.
Immunogen Region	500-550
Specificity	
Immunogen Sequence	