

Anti-IL23A antibody (162-179) {FITC} (STJ501445)

STJ501445

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IL23A (162-179) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.74-0.90 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:250 DB: 1:10, 000 ELISA: 1:10, 000 IP: 0 IHC: 0 ICC: 0 CFM: 0
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51561
Gene Symbol	IL23A
Uniprot ID	IL23A_HUMAN
Immunogen	Synthetic peptide corresponding to C terminal region amino acids 162-179 of Human IL23 P19. Peptide was covalently modified.
Immunogen Region	162-179
Specificity	This antibody does not cross react with other Interleukins or other proteins.
Immunogen Sequence	