

Anti-INHA antibody (1-50) (STJ501461)

STJ501461

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-INHA (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.52-0.75 µg/µl |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 ELISA: 1:20, 000 IP: 1:200 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 3623 |
| Gene Symbol | INHA |
| Uniprot ID | INHA_HUMAN |
| Immunogen | Synthetic peptides corresponding to amino acid region taken within 1-50 and 220-270 on Inhibin alpha gene. |
| Immunogen Region | 1-50 |
| Specificity | |
| Immunogen Sequence | |