

Anti-INHA antibody (315-365) {Biotin} (STJ501468)

STJ501468

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-INHA (315-365) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3623
Gene Symbol	INHA
Uniprot ID	INHA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 315-365 on human alpha Inhibin.
Immunogen Region	315-365
Specificity	
Immunogen Sequence	