

Anti-INSR antibody {Biotin} (STJ501477)

STJ501477

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-INSR is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3643
Gene Symbol	INSR
Uniprot ID	INSR_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	