

Anti-ERN1 antibody (900-977) {Biotin} (STJ501496)

STJ501496

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-ERN1 (900-977) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.59-0.68 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 2081 |
| Gene Symbol | ERN1 |
| Uniprot ID | ERN1_HUMAN |
| Immunogen | Synthetic peptide corresponding to unique amino acid sequence within amino acid region 900-977 on human IRE1A protein. |
| Immunogen Region | 900-977 |
| Specificity | |
| Immunogen Sequence | |