

Anti-F11R antibody (250-299) {Biotin} (STJ501519)

STJ501519

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-F11R (250-299) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Monkey, Mouse, Panda, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52-0.63 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	50848
Gene Symbol	F11R
Uniprot ID	JAM1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-299 on human Junctional Adhesion Molecule A protein.
Immunogen Region	250-299
Specificity	
Immunogen Sequence	