

Anti-KDM4C antibody (500-520) {FITC} (STJ501524)

STJ501524

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KDM4C (500-520) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23081
Gene Symbol	KDM4C
Uniprot ID	KDM4C_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-520 on human JmjC domain-containing histone demethylation protein 3C protein.
Immunogen Region	500-520
Specificity	
Immunogen Sequence	