

Anti-MKI67 antibody (2950-3000) {Biotin} (STJ501532)

STJ501532

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MKI67 (2950-3000) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4288
Gene Symbol	MKI67
Uniprot ID	KI67_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 2950-3000 on human antigen Ki67 isoform 1.
Immunogen Region	2950-3000
Specificity	This antibody will detect both isoform 1 and isoform 2 in human.
Immunogen Sequence	