

Anti-KIF14 antibody (850-960) {FITC} (STJ501539)

STJ501539

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KIF14 (850-960) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.56 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9928
Gene Symbol	KIF14
Uniprot ID	KIF14_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 850-960 on human KIF14 protein.
Immunogen	850-960
Region	
Specificity	
Immunogen	
Sequence	