

Anti-LCT antibody (1080-1130) {FITC} (STJ501566)

STJ501566

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LCT (1080-1130) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.87 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000
	ELISA: 1:10, 000
	IP: 1:200
	IHC: 1:50
	ICC: 1:50
	IF: 1:50
Formulation	CFM: 1:50
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1080-1130 on rat Lactase protein.
Immunogen Region	1080-1130
Specificity	
Immunogen Sequence	