

Anti-LAMP2 antibody (260-310) {FITC} (STJ501572)

STJ501572

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LAMP2 (260-310) is suitable for use in Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	DB, ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.70 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3920
Gene Symbol	LAMP2
Uniprot ID	LAMP2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 260-310 on human lysosomal membrane glycoprotein-2.
Immunogen Region	260-310
Specificity	
Immunogen Sequence	