

Anti-LGR6 antibody (850-920) {Biotin} (STJ501595)

STJ501595

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LGR6 (850-920) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.65 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59352
Gene Symbol	LGR6
Uniprot ID	LGR6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 850-920 on human LRG6 protein.
Immunogen Region	850-920
Specificity	
Immunogen Sequence	