

Anti-LILRB3 antibody (580-631) {FITC} (STJ501608)

STJ501608

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LILRB3 (580-631) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.74 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	102725035/107987462/11025
Gene Symbol	LILRB3
Uniprot ID	LILRB3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 580-631 on human LILRB3 protein.
Immunogen Region	580-631
Specificity	This antibody will detect all three isoforms of human LILRB3 protein.
Immunogen Sequence	