

Anti-LIMK1 antibody (550-647) (STJ501610)

STJ501610

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-LIMK1 (550-647) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgY
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3984
Gene Symbol	LIMK1
Uniprot ID	LIMK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 550-647 on human LIM Kinase 1 protein.
Immunogen Region	550-647
Specificity	
Immunogen Sequence	



Immunohistochemistry of mouse muscle with Anti-LIMK1 primary antibody (STJ501610) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.