

Anti-LOXL1 antibody (570-607) {FITC} (STJ501635)

STJ501635

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LOXL1 (570-607) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.66 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500-1:1,000
Range	ELISA: 1:20,000-1:50,000 IHC: 1:50 IF: 1:50
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16949
Gene Symbol	Lox1
Uniprot ID	LOXL1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 570-607 on rat LOXL1 protein.
Immunogen Region	570-607
Specificity	
Immunogen Sequence	