

Anti-SPACA3 antibody (Full length) (STJ501642)

STJ501642

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SPACA3 (Full length) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	124912
Gene Symbol	SPACA3
Uniprot ID	SACA3_HUMAN
Immunogen	Partially purified full-length lysozyme protein.
Immunogen Region	Full length
Specificity	
Immunogen Sequence	



Immunohistochemistry of human spleen with Anti-SPACA3 primary antibody (STJ501642) at 1:50 dilution in Buffer. Section was treated with Haematoxylin stain. Magnification is 40X.