

## Anti-MAP3K10 antibody (900-950) (STJ501661)

STJ501661

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Chicken polyclonal antibody anti-MAP3K10 (900-950) is suitable for use in Dot Blot, ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	DB, ELISA, IP, WB
<b>Host/Source</b>	Chicken
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.68-0.70 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150
<b>Formulation</b>	
<b>Isotype</b>	IgY
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">4294</a>
<b>Gene Symbol</b>	<a href="#">MAP3K10</a>
<b>Uniprot ID</b>	<a href="#">M3K10_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to the C-terminal taken within amino acid region 900-950 on human MAP3K10 protein.
<b>Immunogen Region</b>	900-950
<b>Specificity</b>	
<b>Immunogen Sequence</b>	