

Anti-MAP3K10 antibody (900-950) {Biotin} (STJ501662)

STJ501662

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MAP3K10 (900-950) is suitable for use in Dot Blot, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	DB, ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.70 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4294
Gene Symbol	MAP3K10
Uniprot ID	M3K10_HUMAN
Immunogen	Synthetic peptide corresponding to the C-terminal taken within amino acid region 900-950 on human MAP3K10 protein.
Immunogen Region	900-950
Specificity	
Immunogen Sequence	