

Anti-MAPK12 antibody (310-360) {FITC} (STJ501689)

STJ501689

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-MAPK12 (310-360) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6300
Gene Symbol	MAPK12
Uniprot ID	MK12_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 310-360 on human Mitogen-activated protein kinase 12.
Immunogen Region	310-360
Specificity	
Immunogen Sequence	