

Anti-MAPK3 antibody (300-350) {FITC} (STJ501704)

STJ501704

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Chicken polyclonal antibody anti-MAPK3 (300-350) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Chicken |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 ELISA: 1:10, 000 IP: 1:150 |
| Formulation | |
| Isotype | IgY |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 5595 |
| Gene Symbol | MAPK3 |
| Uniprot ID | MK03_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 300-350 on human Mitogen-activated protein kinase 3. |
| Immunogen Region | 300-350 |
| Specificity | |
| Immunogen Sequence | |