

Anti-MAPK4 antibody (500-550) {Biotin} (STJ501706)

STJ501706

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-MAPK4 (500-550) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Chicken
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5596
Gene Symbol	MAPK4
Uniprot ID	MK04_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human Mitogen-activated protein kinase 4.
Immunogen Region	500-550
Specificity	
Immunogen Sequence	