

Anti-MAPK7 antibody (700-750) {FITC} (STJ501713)

STJ501713

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MAPK7 (700-750) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5598
Gene Symbol	MAPK7
Uniprot ID	MK07_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-750 on human Mitogen-activated protein kinase 7.
Immunogen Region	700-750
Specificity	
Immunogen Sequence	