

Anti-MAPK8 antibody (370-420) {Biotin} (STJ501715)

STJ501715

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-MAPK8 (370-420) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5599
Gene Symbol	MAPK8
Uniprot ID	MK08_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 370-420 on human Mitogen-activated protein kinase 8.
Immunogen Region	370-420
Specificity	
Immunogen Sequence	