

Anti-MAVS antibody (500-540) (STJ501720)

STJ501720

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MAVS (500-540) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57506
Gene Symbol	MAVS
Uniprot ID	MAVS_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-540 on on human mitochondrial antiviral-signaling protein isoform 1.
Immunogen Region	500-540
Specificity	
Immunogen Sequence	