

Anti-MERTK antibody (900-994) {FITC} (STJ501766)

STJ501766

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MERTK (900-994) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:750-1:2,000 ELISA: 1:10,000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10461
Gene Symbol	MERTK
Uniprot ID	MERTK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-994 on MerTK protein.
Immunogen Region	900-994
Specificity	
Immunogen Sequence	