

Anti-MERTK antibody {Biotin} (STJ501767)

STJ501767

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MERTK is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	10461
Gene Symbol	MERTK
Uniprot ID	MERTK_HUMAN
Immunogen	
Region	
Specificity	
Immunogen Sequence	