

Anti-MTAP antibody (100-200) {Biotin} (STJ501769)

STJ501769

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MTAP (100-200) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution	WB: 1:750
Range	ELISA: 1:10, 000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4507
Gene Symbol	MTAP
Uniprot ID	MTAP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-200 on MerTK protein.
Immunogen Region	100-200
Specificity	
Immunogen Sequence	