

Anti-MYLK2 antibody (200-250) {FITC} (STJ501839)

STJ501839

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MYLK2 (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85366
Gene Symbol	MYLK2
Uniprot ID	MYLK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Myosin light chain kinase 2 protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	