

Anti-MYLK3 antibody (10-60) {FITC} (STJ501842)

STJ501842

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MYLK3 (10-60) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	91807
Gene Symbol	MYLK3
Uniprot ID	MYLK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 10-60 on human Myosin light chain kinase 3 protein.
Immunogen Region	10-60
Specificity	
Immunogen Sequence	