

Anti-MRLN antibody (1-50) (STJ501856)

STJ501856

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MRLN (1-50) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100507027
Gene Symbol	MRLN
Uniprot ID	MLN_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Myoregulin protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-MRLN primary antibody (STJ501856) at 1:50 dilution in buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.