

Anti-PCSK9 antibody (30-60) {FITC} (STJ501875)

STJ501875

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PCSK9 (30-60) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.85-1.25 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100102
Gene Symbol	Pcsk9
Uniprot ID	PCSK9_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 30-60 on the NARC1 protein.
Immunogen Region	30-60
Specificity	
Immunogen Sequence	