

Anti-NEDD4 antibody (840-880) {FITC} (STJ501884)

STJ501884

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NEDD4 (840-880) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:250 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25489
Gene Symbol	Nedd4
Uniprot ID	NEDD4_RAT
Immunogen	Synthetic peptide taken within amino acid region 840-880 on rat Nedd4 protein.
Immunogen Region	840-880
Specificity	
Immunogen Sequence	