

Anti-NKX2-2 antibody (200-250) {Biotin} (STJ501929)

STJ501929

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NKX2-2 (200-250) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.70 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4821
Gene Symbol	NKX2-2
Uniprot ID	NKX22_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 200-250 on human homeobox protein Nkx-2.2.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	