

Anti-NPY2R antibody (300-381) {Biotin} (STJ501964)

STJ501964

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NPY2R (300-381) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64-0.72 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4887
Gene Symbol	NPY2R
Uniprot ID	NPY2R_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-381 on human NPY2R protein.
Immunogen Region	300-381
Specificity	
Immunogen Sequence	