

Anti-NR1H4 antibody (1-50) {Biotin} (STJ501970)

STJ501970

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NR1H4 (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9971
Gene Symbol	NR1H4
Uniprot ID	NR1H4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on huma Bile acid receptor/NR1H4/FXR A3-4 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	