

Anti-D4S234E antibody (140-190) {Biotin} (STJ501973)

STJ501973

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-D4S234E (140-190) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25247
Gene Symbol	Nsg1
Uniprot ID	NSG1_RAT
Immunogen	Synthetic peptide taken within amino acid region 140-190 on human Neuron-specific protein family member 1.
Immunogen Region	140-190
Specificity	
Immunogen Sequence	