

Anti-OCSTAMP antibody (150-200) {FITC} (STJ501995)

STJ501995

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OCSTAMP (150-200) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59-0.63 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10,000 ELISA: 1:10,000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	128506
Gene Symbol	OCSTAMP
Uniprot ID	OCSTP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 150-200 within second intracellular loop of human OC-STAMP protein.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	