

Anti-CDKN2A antibody (Mid-Region) (STJ502021)

STJ502021

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CDKN2A (Mid-Region) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.00-1.25 mg/ml µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1029
Gene Symbol	CDKN2A
Uniprot ID	ARF_HUMAN
Immunogen	Mid-Region epitope synthetic peptide corresponding to unique amino acid sequence on p14ARF protein
Immunogen Region	Mid-Region
Specificity	
Immunogen Sequence	