

Anti-CDKN2A antibody (N-Term) {FITC} (STJ502024)

STJ502024

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CDKN2A (N-Term) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.00-1.25 mg/ml µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1029
Gene Symbol	CDKN2A
Uniprot ID	ARF_HUMAN
Immunogen	N-epitope synthetic peptide corresponding to unique amino acid sequence on p14ARF protein
Immunogen Region	N-Term
Specificity	
Immunogen Sequence	