

## Anti-CDKN2A antibody (Mid-Region) {Biotin} (STJ502025)

STJ502025

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-CDKN2A (Mid-Region) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.00-1.25 mg/ml µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	1029
<b>Gene Symbol</b>	CDKN2A
<b>Uniprot ID</b>	ARF_HUMAN
<b>Immunogen</b>	Mid-Region epitope synthetic peptide corresponding to unique amino acid sequence on p14ARF protein
<b>Immunogen Region</b>	Mid-Region
<b>Specificity</b>	
<b>Immunogen Sequence</b>	