

Anti-CDKN2A antibody (94-111) {Biotin} (STJ502027)

STJ502027

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CDKN2A (94-111) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.00-1.25 mg/ml µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1029
Gene Symbol	CDKN2A
Uniprot ID	ARF_HUMAN
Immunogen	Synthetic peptide corresponding to C-terminal amino acids 94-111 of Human p14ARF.
Immunogen Region	94-111
Specificity	This antibody recognizes p14ARF in human and p19ARF in mouse.
Immunogen Sequence	