

## Anti-CDKN2A antibody (94-111) {FITC} (STJ502028)

STJ502028

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-CDKN2A (94-111) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.00-1.25 mg/ml µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">1029</a>
<b>Gene Symbol</b>	<a href="#">CDKN2A</a>
<b>Uniprot ID</b>	<a href="#">ARF_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to C-terminal amino acids 94-111 of Human p14ARF.
<b>Immunogen Region</b>	94-111
<b>Specificity</b>	This antibody recognizes p14ARF in human and p19ARF in mouse.
<b>Immunogen Sequence</b>	