

Anti-HRAS antibody (100-150) {FITC} (STJ502031)

STJ502031

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HRAS (100-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Chicken, Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.00-1.25 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3265
Gene Symbol	HRAS
Uniprot ID	RASH_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on human GTPase HRAS protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	