

Anti-HRAS antibody (139-189) {Biotin} (STJ502033)

STJ502033

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HRAS (139-189) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Chicken, Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.00-1.25 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10,000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3265
Gene Symbol	HRAS
Uniprot ID	RASH_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 139-189 on human GTPase HRAS protein.
Immunogen Region	139-189
Specificity	
Immunogen Sequence	