

Anti-TP63 antibody (1-100) {FITC} (STJ502049)

STJ502049

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TP63 (1-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59-0.63 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8626
Gene Symbol	TP63
Uniprot ID	P63_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-100 on human p63 protein.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	