

Anti-Adenylate cyclase 1 antibody (C-Term) (STJ502053)

STJ502053

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Adenylate cyclase 1 (C-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.72 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	432530
Gene Symbol	Adcy1
Uniprot ID	ADCY1_MOUSE
Immunogen	Synthetic peptide corresponding to conserved sequence at or near C terminal end of protein that is unique to rat/mouse Adenylyl Cyclase 1 protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	



Western blot with Anti-Adenylate cyclase 1 Antibody (STJ502053) with peptide on at 1:500 dilution in buffer.